

ACTRIS

Aerosol In Situ Data Reporting for 2022

ACTRIS Aerosol In Situ Community Meeting, 25 January 2023



This project receives funding from the European Union's Horizon 2020 research and innovation programme under grant agreements No 871115

ACTRIS Aerosol In Situ Data Reporting for 2022


«The same procedure as last year, Miss Sophie?»

«The same procedure as **EVERY** year, James!»




What to be reported, and how?


- Reporting of all 3 data levels for each instrument: 0, 1, and 2.
- Applies to these ACTRIS aerosol in situ instrument types:
 - Integrating nephelometer
 - Filter absorption photometer (ALL flavours)
 - Mobility particle size spectrometer
 - Condensation particle counter
 - Cloud condensation nucleus counter
 - Aerosol chemical speciation monitor
- Exceptions (only level 2):
 - OC / EC
 - Aerodynamic particle size spectrometer



How are we working?



- National Facilities (NFs / stations) are responsible for data QC!
- Data QC: flagging happens from level 0 (raw as out of instrument) to level 0.a
- Data levels 0.a, 1.a, and 2.a are connected – have to be the same version.
- DC In Situ will create issues in Issue Tracker for each instrument and data level.
- Issues assigned to last year's responsible persons.
- Please notify DC In Situ immediately when responsible person changed!
- Upload finished data submissions to EBAS submission portal.
- **Deadline: 31 May!**
- Helpdesk in late May (use ebas@nilu.no).



Data Quality Control How-To

- Consult [data QC guidelines at ACTRIS ECAC homepage](#).
(separate versions per instrument type)
- Consult templates at [EBAS submission manual](#)
(separate templates per instrument type and data level)
- Pay attention to flag shortlists on template pages
 - Shortlist to make flagging more homogeneous
 - Shortlists vary between instrument type and data level
 - Shortlists slightly longer for level 0 than 1 and 2.
- Use flags 110 and 559 for spike, but instrument worked
 - 110 when you know reason for spike.
 - 559 when you don't know reason for spike.

Group 0: Valid data		
Flag	Validity	Description
000	V	valid data, no flag
Group 1: Exception flags for accepted, irregular data		
Flag	Validity	Description
110	V	Episode data checked and accepted by data originator. Valid measurement
Group 5: Chemical problem		
Flag	Validity	Description
559	V	Unspecified contamination or local influence, but considered valid
Group 6: Mechanical or instrumental problem		
Flag	Validity	Description
640	V	Instrument internal relative humidity above 40%
686	I	Invalid due to zero check. Used for Level 0.
687	I	Invalid due to span check. Used for Level 0.
Group 9: Missing flags		
Flag	Validity	Description
999	M	Missing measurement, unspecified reason

Further Points

- OPSS instruments are NOT on ACTRIS whitelist, but have been reported as ACTRIS in the past
 - Hohenpeissenberg
 - Montseny
 - Pic duMidi
 - OPE
 - La Reunion
 - DEM_Athens
 - Birkenes
- Please continue to report these, but under EMEP and GAW-WDCA
- Please do make sure to have «ACTRIS» in list of frameworks.
- Boundary checks on absorption coefficient data introduced – files overriding check by default will be rejected.
- Normal procedures for AE33: produce and submit all 3 data levels!
- Pay attention to special AE33 metadata elements: filter type (**DOUBLE CHECK!**), nominal flow rate, maximum attenuation, leakage factor zeta, compensation parameters (from setup file)
- Use opportunity: work on AE33 reprocessing in parallel.